

Ultimate Safety, Inc.

4237 Gatecrest, San Antonio, TX 78217

Telephone: (800)248-9241 ~ Fax: (210)599-8250

Ultimate Safety *Fall Prevention*

BedCall

Feature

Benefit

Patented *EdgeSense*

By detecting when the resident positions themselves to the edge of the bed, the controller provides an early warning to staff prior to the resident actually exiting the bed.

Automatic monitoring

The controller automatically begins to monitor the resident once they get into bed eliminating the possibility of someone forgetting to turn on a power switch.

Adjustable time delay

Allows the controller to be set for each resident's unique sleeping habits eliminating false alarms.

Automatic or manual alarm reset

This allows staff to decide whether to require a manual reset of an alarm, or allow the unit to reset if the resident gets back into bed on their own.

Powered by wall adapter

This eliminates high voltage in the controller by using a low voltage adapter that plugs into the wall outlet. No need to worry about batteries running low or when to change them.

Nurse Call Interface

The controller plugs directly into the Nurse Call system and can create an alarm at the central station

Optional local alarm in room

The alarm can be set up to alarm both locally in the room and at the central station

Hold function

The controller can be put on hold by staff so that the resident can be removed from bed without an alarm. Once the resident returns to bed, the unit automatically begins to monitor. The unit remains powered and in standby mode until resident returns.

Volume control

Sound level can be adjusted in room

Ultimate Safety, Inc.

4237 Gatecrest, San Antonio, TX 78217

Telephone: (800)248-9241 ~ Fax: (210)599-8250

Ultimate Safety *Fall Prevention*

ChairCall

Feature

Benefit

Automatic Power-Up

By eliminating the need for an on/off switch, no one can forget to turn the unit on, and the resident can not turn the unit off. This provides the maximum level of security and prevents an accident.

Automatic monitoring

The controller automatically begins to monitor the resident once they sit in the chair eliminating the need to remember to turn on a power switch.

Hold function

The controller can be put on hold by staff so that the resident can be removed from the chair without an alarm. Once the resident returns to the chair, the unit automatically begins to monitor.

Low Battery Alarm

This notifies staff when it is time to change the battery eliminating loss of monitoring.

Nurse Call Input

The controller has a jack for a call cord to plug in that will cause an alarm condition if the resident presses the call button. This allows the resident to summons help while in their chair.

Mounting bracket

The bracket is flexible enough to fit a wide variety of mounting options and can be replaced easily if needed.

Three "AA" batteries

Reduces facilities cost to replace batteries by using AA style as opposed to 9 volt style and provides longer life